promed BODYbelle S

Promed GmbH

Lindenweg 11
D-82490 Farchant
Tel: +49 (0)8821/9621-0
Fax: +49 (0)8821/9621-21
info@promed.de
www.promed.de

RoHS
COMPLIANT

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FI

GETTING STARTED

- A** On/ Off switch and intensity regulation
- B** Mains Plug
- C** Massage head retainer
- D** Massage head
- E** Protective cover cloth

DESCRIPTION OF THE UNIT

Please keep this user manual in a safe place!

Please read the user manual thoroughly before using the appliance for the first time.

Thank you for purchasing the electrical body massager **Promed BODYbelle S**. You have acquired a high quality body massager with many application possibilities. **Promed** is a leading company with decades of experience in personal care, wellness and healthcare.

You can get expert advice for operating the device at all places where you can buy **Promed**-products or you can send us your queries. We may recommend you a responsible consultant. On purchase of a **Promed BODYbelle S**, you get a great device with elegance, luxurious operating comfort in professional quality.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

DESCRIPTION OF THE UNIT

Easy to use, it is ideal for body massage in domestic personal use. Moreover, the attached set with 4 different massage heads, depending on the requirements, helps us in getting rid of excess body fat in key position easily. We, the manufacturer, can in no way be held liable for damage or injury to persons or property arising due to non-observance of this manual.

We wish you every enjoyment with your new **Promed BODYbelle S**. Below we would like to familiarise you with the device. Before first use, please read the direction use.

Scope of delivery:

Device including 4 massage heads and a protective cover cloth.

INFORMATION ON THE UNIT / ACCESSORIES

Please check whether the delivery of the **Promed BODYbelle S** complete and shows no sign of damage. In cases of doubts, do not operate the device and contact your dealer or your service centre. The **Promed BODYbelle S** is provided for an operating voltage of 230V~. Before using the device please verify, whether the voltage required for the device correspond to the local supply (operating voltage 230V~ connected to a socket installed according to the instructions at a supply voltage of 230V~). For the **Promed BODYbelle S**, a high performance motor is used. Adjust the speed of the motor infinitely with the electronic speed control and hence you can specify the massage intensity.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

OPERATING THE UNIT

Before you switch on the device, you must connect the mains plug **(B)** to the mains socket.

The power plug is intended as a separator. Make sure that the power plug is plugged into an easily reachable wall socket.

Operating the device

- A.** Select the desired massage head. Remove the protective cover cloth **(E)** by pulling it from the device. Now, remove the pre-installed massage head **(D)** by pulling it out off the massage head retainer **(C)** . Now insert the desired massage head **(D)** on the massage head retainer **(C)**. Finally, put the protective cover cloth **(E)** back onto the device.
- B.** Now insert the power plug **(B)** into the socket.
- C.** To turn on the device, turn the on/ off switch **(A)** to the left (you must hear a slight click). The controller must be blue now.
- D.** Controlling the massage intensity
To control the massage intensity, you must turn the knob **(A)**.

Guiding to left: Increasing the massage intensity.

Guiding to right: Decreasing the massage intensity.

- E.** After the treatment switch off the massager, pull out the plug and let the device cool.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

SPECIAL QUALITY FEATURES

- Powerful and really deep massage with high intensity, thanks to a very quiet high-performance motor.
- Simple operation, thanks to a variable massage intensity adjustable on a wheel.

In cases of malfunctioning:

- Check that the power supply is correct.
- Check if the on/off switch **(A)** illuminates in blue.
- If the problems have still not been solved after following these steps, please consult your service center immediately.

USAGE

The massage is one of the oldest treatment methods and important part of the physical therapy. The strained muscles are relaxed, blood circulation is improved and the tissue regeneration is encouraged by the mechanical treatment of the outer tissue layers of the body. Thus, the massage has stimulating or soothing effect and can also eliminate neural and muscular pain.

Another effect is the positive effects on the vascular system, internal organs and metabolism. We distinguish between 5 classical massaging techniques: Stroking massage, kneading, tapping and vibration massage as well as friction/slapping. Furthermore, there is another important massage technique such as the manual lymph drainage, reflexology, underwater hydro-massage, brush massage and finger pressure massage (Tai-Chi, Shiatsu, and Acupressure).

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

USAGE

The **Promed BODYbelle S** works on the basis of vibration massage. Here, the massage head is set in vibration by the revolutions of the motor, so that a vibratory motion originates, which is transmitted to body part to be massaged.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

USAGE

- You can massage yourself or a partner. Switch on the device as described above and move it on the body part to be massaged. Do not exert any pressure, but place it gently on the body. The massage is caused solely by the vibration of the device.
- While massaging do not concentrate on the same area (max. 1 minute) for too long, but change the massage area frequently during the use.
- Through the direction of massage you can decide which effect wants to achieve: A massage towards heart relaxes, which stimulates massage away from heart.
- You can also use massage oil for supporting the massage. In this case, the Tai-Chi-massage head is particularly suitable. When using the oils, the protective cover cloth must be cleaned immediately after use.
- An application should not last for more than 10 - 15 minutes.
- Use the device regularly, to configure the positive effect.
- Observe "Safety instructions" while carrying out all massages.
- In principle, the application should always be pleasant. If you feel the pain or find the massage unpleasant, discontinue the use and consult your doctor.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Use of massage heads

	Two sphere-massage head: Acupuncture points and deep massage.
	Four sphere-massage head: Strengthening the blood circulation and the body function.
	Tai-Chi-massage head: Ideal for oil-massage, ensures an enhanced absorption capacity of oil in body.
	Sphere-massage head: Recommended to reduce stress with a neck massage.
	Protective cover cloth: Prevents hair from coming into the machine. Not usable for sphere-massage head.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

MAINTENANCE AND CARE

Cleaning the unit

- Regular maintenance and cleaning ensure optimal results and extends the life of your unit.
- Store the **Promed BODYbelle S** in a safe, cool and dry location when it is not in use.
- To ensure that the safety, reliability and performance of the unit are not adversely affected, enhancements, readjustments, alterations or repairs are not permitted aside from those performed by **Promed** authorized persons.
- Before cleaning the unit, make sure that it is turned off and disconnected from the electricity.
- Use a slightly moist, soft cloth to clean the unit.
- Do not use any aggressive substances such as scouring agents or other detergents as they may corrode the surface.
- Wipe clean with a soft, moist cloth and do not use scouring agents.
- Do not immerse the unit in fluid and keep the unit away from solvents and strong detergents.
- Do not repair the unit. The unit does not have parts that consumers can service.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

Please pay attention to the following before first time operation:

Never place the unit in direct contact with fire, gas, oxygen or hot objects such as hot plates.

Please observe the following during operation:

- A.** Only the supplied massage head must be connected to the device.
- B.** Use the speed controller carefully.
- C.** Switch off the device before changing a massage head.
- D.** You must switch off the device before connecting it to the mains supply (knob **(A)**).
- E.** Power supply to the device: 220V~ - 240V~.
- F.** Take every possible precaution to protect the device from falling or getting damaged otherwise.
- G.** If there are problems in the device, please return it immediately for repair.
- H.** Do not grease or wash the device.
- I.** Only for home use.

- J.** Scrupulous supervision is required when the device is being used in the vicinity of children.

DANGER!

- Do not place the unit in direct contact with water or other fluids.
-  Never place or use the unit in a wet or moist environment.
- Do not touch the unit with wet hands.
- It is not permitted to use the unit outdoors.

WARNING!

- Do not use the device during the pregnancy, on a wound, on inflammations of skin, on skin diseases or in case of cardiovascular diseases.
- Do not use the device within an hour after a meal.
- Do not allow water, oil or other liquids from seeping into device (the massage head can be cleaned by water).
- Make sure that children do not play with the device.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

- Never leave the unit unattended when children or inexperienced persons are in its vicinity.
- Only use the unit for the applications described in this manual.
- Only use attachments and accessories from the manufacturer.
- Never use the unit after the electrical lead has been damaged, when the unit has fallen down or when it has come in contact with water.
- This unit is not designed to be used by people (including children) who have limited physical, sensory, or intellectual abilities or who do not have sufficient experience and/or sufficient knowledge. This does not apply if these people are supervised by a person responsible for their safety or receive instructions from them on how to use this unit.
- Never let children play with the packaging material; there is a risk of suffocation.
- An electrical appliance is not a toy. Children do not recognize the risks that could develop when using electrical appliances. Therefore, operate and store the unit out of the reach of children.
- Do not place or use the unit on hot surfaces or near naked gas flames to ensure that the casing does not melt.
- Never let the unit come into contact with water. Use a moist cloth to clean the unit. Do not use solvents.
- Never carry or pull the unit by the power cord. Use the plug to pull the power cord out of the wall socket. Do not pull on the power cord itself.
- Do not pull the power cord over sharp edges, wedge it in between objects, hang it down or position or hang it over hot objects (e.g., hot plates) or naked flames. Protect the power cord from heat or oil.
- Only use extension cords when you have satisfied yourself beforehand that they are not defect or damaged.
- If the power cable of the appliance gets damaged, have it replaced by an expert immediately.
- Do not store the unit near a wash basin or bath tub as there is the danger that the unit may fall or be pulled into the wash basin or the bath tub.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

- To avoid all risks of injury, blockage or damage to the unit, make sure that the unit does not come into contact with clothing, cords, shoelaces, bristles, wires, threads, hair etc. when in operation.

Note:

The device can heat up slightly. This does not indicate a malfunction. The device can be used without any problem.

DISPOSAL

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

Only for EU countries:

Do not dispose of power tools with the domestic waste!
In accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment and its implementation into national law, nonfunctioning power tools must be collected separately and recycled in an environmentally suitable manner.

This symbol indicates that disposal of this product in the domestic waste is not permitted within the EU. Waste devices contain valuable recyclable material that should be recycled. Recycling also prevents uncontrolled waste disposal from damaging the environment and human health. Therefore, please dispose of all waste devices in the appropriate collection systems or send the device to the place where you bought the device for disposal. They will then recycle the device.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

WARRANTY

This unit has been manufactured with the utmost of care and has been examined in detail before leaving the factory. We therefore provide a warranty of 24 months after the date of purchase when presented with the warranty card issued with this appliance, in accordance with the following conditions.

- Within the warranty period, we will replace all defect parts of the appliance free of charge, including labor costs for repairs under warranty, in cases of verifiable faults in material or manufacture that have occurred with correct usage and have been identified during the warranty period.
 - This warranty is void when the unit has been improperly used (e.g., connection to unsuitable sources of electricity, breakage), the unit has been tampered with (e.g., opening the unit's casing), and when spare parts have been used that have not been approved of by **Promed**. Parts subject to wear are exempt from the warranty. Appliance motors and mobile parts are not included in the warranty.
 - The warranty period begins on the date of purchase. Warranty claims do not affect the warranty duration. Warranty claims must be lodged within the warranty period. Claims arising after cessation of the warranty period cannot be considered.
- The warranty only comes into effect within the scope of these warranty conditions when the date of purchase has been confirmed by the stamp/ signature of the dealer on the warranty card.
 - In cases of warranty or repair claims, please post the complete appliance with the completed warranty card to the customer service point responsible for you. The warranty card can be found at the end of this user manual.
 - We reserve the right to make technical and optical changes as well as changes to the equipment!
 - The following information is vital for problem-free processing:
 1. Original sales slip/receipt or dealer's stamp with the date of purchase
 2. Defect detected
 3. Unit name / type

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Symbol for protection class II appliances

The symbol corresponds to symbol 5172 of the EN 60417-1:2000 Standard. The appliance is double insulated as a protective measure. This prevents conductive parts of the appliance from becoming impaired, something that could lead to malfunctioning of the basic insulation voltage.

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Symbol for WARNING! Indicates sections of the user manual that are relevant for safety!

The symbol corresponds to symbol 0434 of the DIN ISO 7000:2008 Standard. Instructions for an appliance, system or accessory that contain important information for the user that are relevant for safety.

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Pay attention to the symbol for the user manual!

The symbol corresponds to the symbol 1641 in the DIN ISO 7000:2008 Standard. It is imperative that the manual be read before using the unit.

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