beurer

GS 49_BMI

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GB **Glass scale**
Instruction for Use

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Dear customer,

We are delighted that you have chosen to purchase a product from our range. Our name is synonymous with high-calibre and thoroughly tested quality products in the areas of thermals, weight, blood pressure, body temperature, pulse, gentle therapy + beauty, massage and air.

Please read through this instruction manual carefully and retain it for future reference; make it available to other users and observe the information.

Best wishes,
Your Beurer team

1. Good to know – How to weigh correctly

As far as possible, weigh yourself at the same time everyday (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.

In order to calculate your height simply and accurately, the scale uses a height calculator, which measures your height. It has a measurement range of between 50 cm and 2 m.

! 2. Important information – retain for future reference!

Warning – for your safety

- Warning, do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!
- Protect the scale from knocks, damp, dust, chemicals, drastic changes in temperature and nearby sources of heat (ovens, heaters).
- Repairs may only be carried out by Beurer Customer Services or authorised retailers. Before each repair query, first check the batteries and replace them if necessary.
- Each scale is made in accordance with EC Directive 2004/108 + amendments. If you should have further questions about using our devices, please contact either your retailer or Customer Services.
- The device should be cleaned with a damp cloth from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Storage: Do not place any objects on the scale when they are not being used.
- The device is suitable only for self-testing, not for medical or commercial purposes.
- When supplied to the customer, the scale is set to weigh in “kg”. To the rear of the scale, there is a button which you can use to change to “pounds” and “stones”.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile telephones).

3. Use

This scale is used to calculate people's weight and BMI. For this purpose, this scale has a weighing range of up to 150 kg. If the scale is loaded with a higher weight (>150 kg), “-----” will appear in the display. Your height is necessary for calculation of your BMI. If you simply wish to calculate your height, please use the provided height calculator. If your height is outside of the measurement range (50 cm – 2 m), “Err1” will appear. If calibration is incorrect, “Err2” will appear.


3.1 Operation

If there is one, remove the plastic insulating strip by pulling it from the battery compartment lid of the scale and the height calculator and remove the protective film from the batteries, then insert them the correct way round. If the scale does not function, remove the batteries completely and then replace them.

Your scale is equipped with a “low battery indicator”. If you operate the scale with flat batteries, “Lo 2” will appear on the display and the scale will automatically switch itself off. In this case, the batteries in the scale must be replaced (2 x CR2032). If “Lo 1” appears in the display, the batteries of the height calculator must be changed. (2 x AAA 1.5 V micro).


If the scale does not function, change the batteries. Change the height calculator batteries if it does not beep when taking height measurements (exception: the height calculator does not beep when you hold it horizontally).

The empty, completely flat batteries should be disposed of by means of specially designated collection boxes, recycling points or via electronics retailers. You are legally required to dispose of the batteries.

i **Tip:** You will find these symbols on batteries containing contaminants: Pb = batteries contain lead, Cd = batteries contain cadmium, Hg = batteries contain mercury.

Please dispose of the device in accordance with the EC directive on old electric and electronic devices 2002/96/EC – WEEE (“Waste Electrical and Electronic Equipment”). If you have queries, please contact the appropriate local authorities.

3.2. Only measuring weight

Place the scale on an even, hard surface (no carpet); a hard surface is crucial to accurate measurements.





Step carefully onto the scale, distributing your weight equally across both legs. The scale begins to take measurements straight away. After the display has flashed twice, the final measurement is displayed (Fig. 1). If you step off the weighing area, the scale switches off after 10 seconds.

Fig. 1

3.3 Calculating weight and BMI

- Stand still on the scale with equal weight distribution until your weight is displayed (Fig. 2). The scale flashes twice to indicate that your weight has been measured.
- Position the height calculator on the middle of your head, with the thin end pointing forwards and the button pointing upwards.
- Press the button on the height calculator once (do not keep the button pressed down). It will then begin to beep (Fig. 3).
- Look straight ahead.
- Use the height calculator to find a horizontal position. Keep the thicker end on your head (Fig. 4).
- As soon as you have found a horizontal position, the beeping will stop. Remain in this position for about 5 seconds. (Fig. 5)
- The scale beeps and signals that the height has been measured (Fig. 6). The measurement has now been taken and you may put the pole aside.
- The scale alternately displays your weight, your BMI and your height.
Tip: The BMI can only be calculated if the weight has been correctly measured.
- The arrow functions on the right-hand edge of the display and the colour chart on the device show which area the recorded BMI values fall into.
- You can now step off the scale. The scale switches off automatically after 10 seconds.

i Tip:

- The height calculator stops beeping automatically after about 30 seconds.
- Also, the beep does not sound as long as the height calculator is in a horizontal position.
- The height measurement is temperature dependent. In order to achieve precise results, use the scale and height calculator at room temperatures between 20° C and 30° C.
- When using the scale for the first time or after moving the scale, touch the weighing area first without standing on the scale. After turning off the display, the normal measurement can be taken.

4. Analysing the results

The body mass index (BMI) is a popular method used to interpret body weight. It analyses your weight in relation to your height. This is the formula:

$$\text{BMI} = \frac{\text{weight (kg)}}{(\text{height (m)})^2}$$

The WHO (World Health Organisation) has classified the values for adults over 20 years old as follows:

Colour chart	Category	BMI (KG/M ²)	
blue	Severely underweight	< 16.00	< 18.50 means underweight
blue	Underweight	16.00 - 19.99	
blue	Slightly underweight	17.00 - 18.49	
green	Normal weight	18.50 - 24.99	
Orange	Overweight	25.00 - 29.99	> 25 means overweight
Red	Class I obesity	30.00 - 34.99	> 30 means obese
Red	Class II obesity	35.00 - 39.99	
Red	Class III obesity	> 40.00	

The interpretation which appears on the display is valid for adults over 20 years old. There are separate tables for children.

To avoid obesity and any accompanying illnesses, you should aim to stabilise or slightly reduce your weight if your BMI is between 25 and 29.9.

However, if your values lie in the “overweight” category you are recommended to steadily lose weight if:

- you suffer from existing weight-related illnesses such as high blood pressure or type 2 diabetes
- you have uneven fat distribution (gaining weight around the stomach – the “apple figure”)
- you suffer from existing diseases that are exacerbated by being overweight (such as back or joint problems or asthma)
- you are exposed to great psychosocial pressure

i **Tip:** For children and young people of 5-19 years old, the following guideline values are classified by the WHO:

Signs for girls:

Age	Severely underweight	Underweight	Normal	Overweight	Obese
5	< 11.7	11.7 – 12.7	12.8 – 17.0	17.1 – 19.2	< 19.2
6	< 11.7	11.7 – 12.7	12.8 – 17.3	17.4 – 19.7	> 19.7
7	< 11.9	11.9 – 12.9	13.0 – 17.7	17.8 – 20.5	> 20.5
8	< 12.1	12.1 – 13.1	13.2 – 18.3	18.4 – 21.4	> 21.4
9	< 12.4	12.4 – 13.4	13.5 – 19.0	19.1 – 22.5	> 22.5
10	< 12.7	12.7 – 13.8	13.9 – 19.8	19.9 – 23.6	> 23.6
11	< 13.1	13.1 -14.3	14.4 – 20.7	20.8 – 24.9	> 24.9
12	< 13.6	13.6 – 14.9	15.0 – 21.7	21.8 – 26.1	> 26.1
13	< 14.0	14.0 – 15.4	15.5 – 22.7	22.8 – 27.2	> 27.2
14	< 14.3	14.3 – 15.8	15.9 – 23.5	23.6 – 28.2	> 28.2
15	< 14.6	14.6 – 16.2	16.3 – 24.1	24.2 – 28.8	> 28.8
16	< 14.7	14.7 – 16.3	16.3 – 25.5	25.6 – 29.3	> 29.3
17	< 14.7	14.7 – 16.4	16.5 – 24.8	24.9 – 29.5	> 29.5
18 - 19	< 14.7	14.7 – 16.5	16.6 – 25.0	25.1 – 29.7	> 29.7

Signs for boys:

Age	Severely underweight	Underweight	Normal	Overweight	Obese
5	< 12.1	12.1 – 13.0	13.1 – 16.7	16.8 – 18.5	> 18.5
6	< 12.2	12.2 – 13.1	13.2 – 17.0	17.1 – 19.0	> 19.0
7	< 12.3	12.3 – 13.3	13.4 – 17.4	17.5 – 19.6	> 19.6
8	< 12.5	12.5 – 13.5	13.6 – 17.9	18.0 – 20.4	> 20.4
9	< 12.8	12.8 – 13.7	13.8 – 18.4	18.5 – 21.3	> 21.3
10	< 13.0	13.0 – 14.0	14.1 – 19.1	19.2 – 22.4	> 22.4
11	< 13.4	13.4 – 14.4	14.5 – 19.9	20.0 – 23.5	> 23.5
12	< 13.8	13.8 – 14.9	15.0 – 20.8	20.9 – 24.7	> 24.7
13	< 14.2	14.2 – 15.4	15.5 – 21.7	21.8 – 25.8	> 25.8
14	< 14.7	14.7 – 16.0	16.1 – 22.6	22.7 – 26.9	> 26.9
15	< 15.1	15.1 – 16.5	16.6 – 23.5	23.6 – 27.8	> 27.8
16	< 15.4	15.4 – 16.9	16.9 – 24.2	24.2 – 28.6	> 28.6
17	< 15.7	15.7 – 17.3	17.4 – 24.9	25.0 – 29.2	> 29.2
18	< 15.8	15.8 – 17.5	17.6 – 25.4	25.5 – 29.7	> 29.7
19	< 15.9	15.9 – 17.6	17.7 – 25.4	25.5 – 29.7	> 29.7

Not all those classified as overweight are automatically at greater risk of illness. If you have any doubts, consult your doctor.

For those classified as overweight, further criteria such as waist measurement should also be drawn upon. For it is not only the BMI which serves as a guideline; fat distribution also plays a decisive role. Investigations have shown that the health risk is increased by excessive fat distribution to the abdomen. The guideline value here is waist circumference, taken from between the bottom rib bone and the top of the pelvis.

Men are at higher risk at waist measurements >94 cm, women at >80 cm.

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