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GETTING STARTED	p. 19	NORMALIZED SYMBOLS	p. 181/182
CONTENTS	p. 20	CE DIRECTIVES	p. 185/186
SYSTEM COMPONENTS	p. 21	TECHNICAL SPECIFICATIONS	p. 187/188
DESCRIPTION OF THE UNIT	p. 22	TRANSPORT/STORAGE/OPERATING CONDITIONS	p. 189/190
INTRODUCTION	p. 23	ACCESSORIES	p. 191/192
INDICATIONS / CONTRAINDICATIONS	p. 24	CLINICAL RECOMANDATIONS	p. 195/196
WARNINGS / PRECAUTIONARY MEASURES	p. 25/26	ANNEX A PLACEMENT OF ELECTRODES	p. 197
OPERATING THE UNIT	p. 27-30	WARRANTY CARD	p. 198
PROGRAM	p. 31/32	MANUFACTURER + FURTHER INFORMATION	p. 0
SAFETY INSTRUCTIONS	p. 33		
TROUBLESHOOTING	p. 34		
CARE, MAINTENANCE AND STORAGE	p. 35		
DISPOSAL	p. 35		
WARRANTY & SERVICING	p. 36		

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Your **Promed T-6 - package** contains the following equipment:

- |                             |       |
|-----------------------------|-------|
| 1. T-6 Unit                 | 1 PCS |
| 2. Cable                    | 2 PCS |
| 3. Electrode pads (40*40mm) | 4 PCS |
| 4. User manual              | 1 PCS |
| 5. AAA battery              | 4 PCS |
| 6. Carrying case            | 1 PCS |
| 7. Belt Case                | 1 PCS |

Please read the user manual thoroughly before using the device for the first time.

Thank you for purchasing the **Promed T-6**. You have purchased a high-quality product designed for personal hygiene and wellbeing purposes. **Promed** is a leading company with decades of experience in personal hygiene, wellness and health areas.

The **Promed T-6** unit has been designed and manufactured in accordance with the Medical Device Directives 93/42/EEC and 2007/247/ EEC to guarantee quality during use. The device may be used once this user manual has been read.

We, as the manufacturer, cannot be made liable in any way for injury or damages to people or objects that arise from failure to comply with this user manual. We wish you a lot of enjoyment with your new **Promed T-6**. We would like to familiarize you with the TENS unit in the following sections. Please read the user manual thoroughly before using the device for the first time.

You can receive professional advice wherever **Promed** products are sold or get in touch with us if you have any questions. We can provide you with the name of the representative responsible for you.

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Further informations can be found on our homepage:

**[www.promed.de](http://www.promed.de)**

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**Promed T-6** is a two output channels TENS device and highly effective in relieving pain. It is intended for symptomatic relief and management of chronic (long term) intractable pain and or as an adjunctive treatment in the management of post surgical and post traumatic acute pain problem. TENS stands for Transcutaneous Electrical Nerve Stimulation is recognized as a safe and highly effective method of relieving pain. TENS is now regularly recommended by doctors, physiotherapists and pharmacists throughout the world. Before using, please read all the instructions in this user manual carefully and keep it safe for future reference.

### What is pain?

Pain is the body's own early warning system. Feeling pain is important, as it indicates an abnormal condition within the body and serves as a warning to us before any further damage or injury occurs. However, long-term, persistent pains – often known as chronic pain – serve no obvious useful purpose once a diagnosis has been made. TENS was developed to soothe or eliminate certain types of chronic and acute pain.

### We differentiate between two types of pain:

#### • Acute pain

as a chief symptom can often help the doctor with diagnosis and the acute pain has a protective function for the patient.

#### • Chronic pain

can often become part of the illness in itself. A patient suffering from chronic pain will often suffer for years and experiences changes in his/her personality structure.

### EXPLANATION OF TENS

Transcutaneous Electrical Nerve Stimulation (TENS) is a non-invasive, drug free method of controlling pain. TENS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. TENS does not cure any physiological problem; it only helps control the pain. TENS does not work for everyone; however, in most patients it is effective in reducing or eliminating the pain, allowing for a return to normal activity.

### How does TENS work?

**Promed T-6** works by passing harmless electrical signals into the body from the electrode pads. This relieves pain in two ways:

- Firstly, it blocks the body's pain signals which are normally transmitted from the area of damage through the nerve fibres to the brain - TENS interrupts these pain signals.
- Secondly, TENS stimulates the body's production of endorphins – its own natural painkillers.

### Possible TENS applications

Generally the **Promed T-6** shall be used for the following medical indications or the following complaints as a TENS unit for treatment:

- For symptomatic relief of chronic intractable pain
- For post traumatic pain (acute pain occurring)
- For post surgical pain (pain caused by surgery)

*Be sure to consult a doctor before using the unit!*

*In the USA, the sale of this device is restricted by law to doctors, or subject to a doctor's prescription.*

Please carefully read and understand the following warnings and cautions to ensure the safe and correct use of this device and to prevent injury.

### Indications

Transcutaneous Electrical Nerve Stimulation (TENS) can be prescribed by doctors for the symptomatic reduction and treatment of chronic (long-term) pain and to treat postoperative or post-traumatic pains.

### Contraindications

Do not use this device without consulting a doctor if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.

If you apply to one of the following cases always be sure to consult a doctor before you apply the **Promed T-6** and clarify the use of the device with him.

- In case of persistent pain symptoms, despite therapy
- Taking a strong analgesic or local anesthetic
- In Infectious Diseases
- By circulatory disorder (thrombosis and embolism)
- For sensory disturbances (numbness)
- When used in infants, toddlers and children

- In pregnancy
- With Psychosis
- If bleeding tendency
- With Cancer
- At extreme sensitivity to or fear of electricity
- In patients with metallic implants
- For heart problems, especially cardiac arrhythmias
- Prior to each electrode placement, that apply current to the carotid sinus region (anterior neck)
- Prior to each electrode placement, that apply current transcebrally (through the head)
- If undiagnosed pain symptoms
- For treatment of the eyelids
- Serious arterial circulatory problems (embolism) in the lower limbs
- The device should not be used for symptomatic local pain relief unless etiology is established or unless a pain syndrome has been diagnosed
- When cancerous lesions are present in the treatment area
- At swollen, infected, inflamed areas or skin eruption (e.g. phlebitis, thrombophlebitis, varicose veins, etc.)
- On-existence of a demand-type cardiac pacemaker or any implanted defibrillator
- In areas of the body with poorly enervated areas
- With epilepsy
- For abdominal or inguinal hernia

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**⚠ Warnings**

- TENS devices must be kept out of the reach of children.
- The safety of using TENS devices during pregnancy or delivery has not been proven.
- TENS is not effective in treating pain of the central nervous system (headaches).
- Should TENS treatment fail to work or cause discomfort, stimulation should be discontinued until a doctor has reassessed the situation.
- Always switch the TENS unit off before attaching or removing the electrodes.
- Never position the electrodes over the eyes, in the mouth or within the body.
- TENS devices do not possess any healing properties.
- TENS is a method for treating symptoms and as such only suppresses the perception of pains that would otherwise serve as a defence and warning mechanism for the body.
- Do not apply stimulation while driving, operating machinery, or during any activity in which electrical stimulation can put you at risk of injury.
- If you are in the care of a physician, consult with your physician before using this device.
- If you have had medical or physical treatment for your pain, consult with your physician before using this device.
- If your pain does not improve, becomes more than mild, or continues for more than five days, stop using the device and consult with your physician.
- Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not apply stimulation when in the bath or shower.
- Do not apply stimulation while sleeping.
- Do not use the device on children, if it has not been evaluated for pediatric use.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Apply stimulation only to normal, intact, clean, healthy skin.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

### Precautionary measures / Negative reactions Precautions

- TENS is not a substitute for pain medications and other pain management therapies.
- TENS devices have no curative value.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.
- The long-term effects of electrical stimulation are unknown.
- You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium.
- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.
- Consult with your physician prior to using the device after a recent surgical procedure, because stimulation may disrupt the healing process.
- Use caution if stimulation is applied over the menstruating or pregnant uterus.
- Use caution if stimulation is applied over areas of skin that lack normal sensation.
- Use this device only with the leads, electrodes, and accessories recommended by the manufacturer.

### Negative Reactions

- Long-term use can occasionally lead to skin irritation in the area where the electrodes were placed.
- The effectiveness of the treatment depends greatly upon the patient being treated by someone who is qualified in dealing with patients suffering from pain.
- Potential negative reactions include skin irritation and burns caused by the electrodes.
- You may experience skin irritation and burns beneath the stimulation electrodes applied to your skin.
- You may experience headache and other painful sensations during or following the application of electrical stimulation near your eyes and to your head and face.
- You should stop using the device and should consult with your physician if you experience adverse reactions from the device.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

**Before you start:**

- Make sure the batteries are installed correctly.
- Connect the cable to the electrode pads and then connect to the unit.
- Attach the electrode pads to the area of pain.
- Do not use the electrode pads if it is scratched or damaged in any way.

**Check/Replace the battery**

1. Remove the Belt Clip (**J**) by sliding it down the unit.
2. Open the battery compartment cover (**I**).
3. Insert 4 batteries (type: AAA) into the battery compartment. Make sure you are installing the batteries correctly. Be sure to match the positive and negative ends of the batteries to the markings in the battery compartment of the device.
4. Recover the battery compartment (**I**) cover if you want to use it.
5. Bring on the Belt Clip (**J**) back by sliding it upwards the unit.

**Caution:**

1. Remove the batteries if the device is not in use for long periods of time.
2. Do not mix old and new batteries or different types of batteries.


3. Warning: If batteries leak and come into contact with the skin or eyes, wash immediately with copious amounts of water.
4. Batteries must be handled by an adult. Keep batteries out of the reach of children.
5. Only batteries of the same or equivalent type are recommended.
6. Do not use rechargeable batteries.
7. Remove exhausted batteries from the unit.
8. Dispose of batteries safely according to battery manufacturer's instructions.

**Connect the electrode pads to the cables:**

Hold the cable plug and insert into electrode pad connections. Make sure there are no bare metal of the pins exposed.

**Caution:**

Always use the cable that supplied with manufacturer or distributor and use the electrode pads with CE mark, or are legally marketed in the US under 510 (k) procedure.

**Connect cable to unit**

Before proceeding to this step, ensure the device is completely switched off, Hold the cable plug and insert into the socket (**B**) on the top of the unit.

**Caution:**

Do not insert the plug of the patient lead wire into any AC power supply socket.

**Place electrode pads on skin**

Apply electrode pads to the area where the pain is felt (refer to Annex A: placement of electrodes). Before applying electrodes, be sure the skin surface over which electrode pads are placed is thoroughly clean and dry. Make sure the electrode pads are pressed firmly to the skin and make good contact between the skin and the electrode pads. Place the electrode pads over the skin, attach them properly, firmly, and evenly.


**Caution:**

1. Before applying the electrode pads to the body, ensure the skin surface is clean and free from lotions or moisturizers.
2. Do not switch the device on when self-adhesive electrode pads are not positioned on the body.
3. Never remove the pads from the skin whilst the device is still switched on.
4. For hygiene reasons we recommend that the electrode pads are replaced every 30 days.
5. Please make sure the device is completely turned off before connecting.
6. It is recommended that, at minimum, 40mm x 40mm self-adhering based, square electrodes are used at the treatment area.
7. For hygienic reasons, each patient should use their own set of electrodes.

**Turn the device on:**


Push up the “Power ON/OFF” switch (A) on the right side of unit to turn the device on, the following screen will appear:

**Select the program:**

There are 11 (P01~P11) preset programs and 4 (P1~P4) user programs available for you to select from. The details of the programs, please refer to section programs.

Press the “MODE” button (E), you can select the program you need. And the LCD will display like following:

### Set the parameter

In the user programs (P1~P4), the user can set the pulse rate, pulse width and treatment time by pressing the "SET" button (F). The LCD will be display the current parameter which you want to set like following:

To set the parameter you can press "▼" (G) or "▲" (H) button.

### Remark:

The set range of pulse rate, pulse width and treatment time in the user program, please refer to section programs.

### Adjust the intensity and Start to treatment:

There are two intensity adjustment knobs (C) for each channel on the side of the unit. By turning the controls up, the output intensity will be increased, when the output intensity above 1mA, the device will be started to work. The maximum output intensity is 100mA.


### Remark:

The **Promed T-6** has special backlight function:

- In the normal state the backlight is green;
- If the output intensity is over 65mA, the backlight will be changed blue and the user should be taken care.
- If the electrode pad connection in wrong way or not connect the electrodes, when the output intensity over 10mA, the backlight will change to red and flash, and the intensity will be reseted automatically.

### Lock function

In order to use the device safely, after you adjusted the output intensity, please press the "SET" button (F) and "▲" button (H) at the same time, the lock function mark "🔒" (9) will show up on the LCD (The device will be locked). This is a safety feature to prevent accidental changes to your settings and to prevent accidental increases the intensity level. If you want to unlock, Press the "SET" button (F) and "▼" button (G) at the same time, the lock function mark "🔒" (9) will disappear.

### Check the recording memory

1. The **Promed T-6** device can record 30 times treatment data. In the waiting or setting state, to check the memory data, press the "SET" button (**F**) and hold on for 5 seconds, the device will enter the memory state. First, the last therapeutic program recoding will display on the LCD, you can press the "SET" button (**F**) again to check the treatment parameter of this program (treatment time, pulse rate and pulse width).
2. Press "▼" (**G**) or "▲" (**H**) button to check other therapeutic program memory.
3. If you want to back the waiting or setting state, press the "MODE" button (**E**) or wait for 30 seconds without any operation in the panel.
4. To clear the memory, press the "SET" button (**F**) and hold on for 5 seconds, the symbol of "D" will be flash on the LCD to remind the user to delete the memory. You can press the "SET" button (**F**) again to delete the memory.

### Low battery indicator:

When the low power indicator (**B**) flashes, the batteries should be replaced with new ones as soon as possible. However, the device will continue to operate for several more hours. When the device power off, the current state will be saved automatically.

### Turn the device off

Push "Power ON/OFF" switch (**A**) down to turn off the device.

### Caution

If there is no operation in the panel in waiting state, you will hear a "DI" sound after one minute and hear two "DI" sounds after two minutes; and then the device will automatically turn into the power saving state after three minutes, LCD will close the display and you will hear three times "DI" sounds. You can press any button in the panel to activate the display.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

There are 11 (P01~P11) preset programs and 4 (P1~P4) user programs available for you to select from. The Treatment time, pulse rate and pulse width can be set in user program (P1~P04), the setting range of the parameter please refer to below. And the electrode pad positions please see Annex A Placement of electrodes.

**preset programs**

Program	Pulse rate	Pulse width	Waveform	Treatment time	Indication for use / Electrode placement no.
P01	80 Hz	180 $\mu$ S	Con.	30 Min	Neck pain 2 Shoulder pain 3 Elbow pain 7 Rheumatic pain 9 Lumbago 11,12 Menstrual pain 13 Phantom limb pain 14 Hip pain 16 Osteoarthritis pain in the knee 18,19 wound healing 20
P02	Burst rate: 2 Hz Frequency fixed: 100 Hz	180 $\mu$ S	Burst	30 Min	Cervical rhizopathy 4 Central pain 5,6 Sciatica 15 Knee pain 19
P03	80 Hz	50 $\mu$ S -180 $\mu$ S	Pulse width modulation	30 Min	Trapezius pain 10, 21 Lumbago 11,12
P04	2 Hz / 15 Hz	180 $\mu$ S	Alternating	30 Min	Neck pain 2 Shoulder pain 3
P05	2 Hz / 80 Hz	180 $\mu$ S	Alternating	30 Min	Lumbago 11,12 Hip pain 16
P06	2 Hz / 100 Hz	200 $\mu$ S /150 $\mu$ S	Alternating	30 Min	Thigh pain 17

P07	10 Hz	180 $\mu$ S	Continuous	30 Min	Postoperative or Chemotherapy-induced Nausea	8
P08	80 Hz	60 $\mu$ S	Continuous	30 Min	Facial pain Neck pain	1 2
P09	2 Hz	250 $\mu$ S	Continuous	30 Min	See program P07	
P10	2 Hz - 80 Hz	200 $\mu$ S - 100 $\mu$ S	Modulation	30 Min	Shoulder pain Trapezius pain Lumbago Thigh pain	3 10, 21 12 22
P11	100 Hz	75 $\mu$ S	Continuous	30 Min	See program P02	

**user programs**

<b>Program</b>	<b>Pulse rate</b>	<b>Pulse width</b>	<b>Waveform</b>	<b>Treatment time</b>	<b>Indication for use / Electrode placement no.</b>
[P1]	Burst rate: 0.5 Hz - 5 Hz Frequency fixed: 100 Hz	30 $\mu$ S - 300 $\mu$ S	Burst	1 Min - 60 Min /Continue	See program P02
[P2]	1 Hz - 150 Hz	30 $\mu$ S - 300 $\mu$ S	Continuous	1 Min - 60 Min /Continue	See program P01
[P3]	1 Hz - 150 Hz	30 $\mu$ S - 300 $\mu$ S	Pulse width Modulation	1 Min - 60 Min /Continue	See program P03
[P4]	1 Hz - 150 Hz	30 $\mu$ S - 300 $\mu$ S	Pulse rate mo- dulation	1 Min - 60 Min /Continue	See program P03

## SAFETY INSTRUCTIONS

### Pay attention to the following before first time operation:

- Only for home use.
- Scrupulous supervision is required when the device is being used in the vicinity of children.
- Never place or use the unit in a wet or moist environment.
-  Do not use the unit under water, e.g., under the shower.
- Never place the unit in direct contact with fire, gas or oxygen, or hot objects such as hot plates.
- Take every precaution possible to ensure that the unit does not fall down or get damaged in any other way.
- If problems do arise with the unit, take it immediately to get it repaired.
- Do not lubricate or wash the unit.

### **Danger!**

- Do not place the unit in direct contact with water or other fluids.
- It is not permitted to use the unit outdoors.
- Never touch the unit with wet hands.
- Do not store the unit near a wash basin or bath tub as there is the danger that the unit may fall or be pulled into the wash basin or the bath tub.

### **Warning!**

- Never leave the unit unattended when children or inexperienced persons are in the vicinity of the unit.
- Make sure that children do not play with the device.
- Only use the unit for the applications described in this manual.
- This unit is not designed for use by people (including children) who have limited physical, sensory, or intellectual abilities or who do not have sufficient experience and/or sufficient knowledge. This does not apply if these people are supervised by a person responsible for their safety or receive instructions from them on how to use the unit.
- Do not use the unit underneath blankets or cushions as excessive heat may cause fire, injury or an electrical shock.
- Do not carry the unit by its cable and do not use the cable as a handle.
- It is not permitted to use the unit outdoors.
- When you have finished treatment, place the intensity controllers **(2)** in the "OFF" position.
- Never let children play with the packaging material; there is a risk of suffocation.

DE

EN

FR

IT

ES

NL

RU

PL

SE

FI

Problem	Possible Cause	Solution
Displays fail to light up	Battery contact failure	1. Try fresh batteries.
		2. Ensure batteries are inserted correctly. Check the following contacts: <ul style="list-style-type: none"> <li>• All contacts are in place.</li> <li>• All contacts are not broken.</li> </ul>
Stimulation weak or cannot feel any stimulation	Electrodes 1. Dried out or contaminated 2. Placement	Replace and re-connect
	Lead wires Old/worn /damaged	Replace
	The intensity is too weak.	Use a higher intensity level.
Stimulation is uncomfortable.	Intensity is too high	Decrease intensity.
	Electrodes are too close together	Reposition the electrodes.
	Damaged or worn electrodes or lead wires	Replace
	Electrode active area size is too small.	Replace electrodes with ones that have an active area no less than 16.0 cm <sup>2</sup> (4cm * 4cm).
	Mayn't operate the device according to the manual.	Please check the manual before use.
Intermittent output	Lead wires	1. Verify connection is secure. Insure firmly.
		2. Turn down the intensity. Rotate lead wires in socket 90°. If still intermittent, replace lead wire.
		3. If still intermittent after replacing lead wire, a component may have failed. Call the repair department.
		4. Some programs will seem intermittent. This is expected.

Stimulation is ineffective.	Improper electrode and applicator placement Unknown.	Reposition electrode and applicator-Contact clinician.
The skin becomes red and/or you feel a stabbing pain	Use the electrodes on the same site every time.	Re-position the electrodes. If at any time you feel pain or discomfort stop use immediately.
	The electrodes aren't stuck onto the skin properly.	Ensure the electrode is stuck securely on the skin.
	The electrodes are dirty.	Clean the electrode pads with a damp, lint free cloth or replace new electrode pads.
Output current stops during therapy	The surface of the electrode was scratched.	Replace new electrode.
	The electrode pads come off the skin.	Turn off the device and stick the electrode pad firmly to the skin.
	The cable is disconnected.	Turn off the device and connect the cable
	The power of the batteries has been exhausted.	Please replace them with new batteries.

- Do not place the unit in direct contact with water or other fluids.
- Switch the unit off when it is not in use.
- Before putting the **Promed T-6** into storage for a longer period of time, remove the batteries from the compartment. Leaking batteries can damage the unit.
- Keep the unit and its accessories in a cool, dry place in the case provided.
- Handle the electrodes by their edges when removing them. To avoid damaging them, never pull at the electrode leads themselves.
- Do not make any sharp kinks in the connecting leads or electrodes.
- After use, stick the electrode pad onto the protective plastic film.
- Do not expose the device to direct sunlight and protect it against dirt and moisture.
- Never place any heavy objects on the device.
- You can clean your **Promed T-6** by wiping it over carefully with a cloth moistened with mild soapy water. You may also use isopropyl alcohol or a soap solution. Household detergents and cleaning products are not suitable.
- The unit must not be returned to your dealer or the manufacturer for inspection and recalibration. If you need such inspections or recalibrations because of your internal Quality Management System don't hesitate to contact Promed.
- If you should experience other problems, consult your dealer, returning the unit if necessary. Never attempt to repair a fault yourself.

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

**Only for EU countries:** Do not throw power tools in with the domestic waste! In accordance with the European Directive 2002/96/EEC on waste electrical and electronic equipment and its implementation into national law, nonfunctioning  power tools must be collected separately and recycled in an environmentally suitable manner.

This symbol indicates that disposal of this product in the domestic waste is not permitted within the EU. Waste devices contain valuable recyclable material that should be recycled. Recycling also prevents uncontrolled waste disposal from damaging the environment and human health. Therefore, please dispose of all waste devices using the appropriate collection systems or send the device to the place of purchase for disposal. The place of purchase will then recycle the device.

#### **Rechargeable batteries/batteries:**

Do not throw your rechargeable batteries/batteries in the domestic waste, in fire, or in water. Rechargeable batteries/batteries should be collected, recycled, or disposed of in an environmentally friendly manner.

#### **Only for EU countries:**

According to the directive 91/157/EEC, defect or used rechargeable batteries/batteries must be recycled. Waste rechargeable batteries/batteries can be directly handed to:  
Promed GmbH, Lindenweg 11, D-82490 Farchant.

This unit has been manufactured with the utmost of care and has been examined in detail before leaving the factory. We therefore provide a warranty of 24 months after the date of purchase when presented with the warranty card issued with this appliance and in accordance with the following conditions:

- Within the warranty period, we will replace all defect parts of the unit free of charge, including labor costs for repairs under warranty, in cases of verifiable faults in material or manufacture that have occurred with correct usage and have been identified during the warranty period.
- This warranty is void when the unit has been improperly used (e.g., connection to unsuitable sources of electricity, breakage), the unit has been tampered with (e.g., opening the unit's casing), and when spare parts have been used that have not been approved of by **Promed**. Parts subject to wear are exempt from the warranty. Appliance motors and mobile parts are not included in the warranty.
- The warranty period begins on the date of purchase. Warranty claims do not affect the warranty duration. Warranty claims must be lodged within the warranty period. Claims arising after cessation of the warranty period cannot be considered.
- The warranty only comes into effect within the scope of these warranty conditions when the date of purchase has been confirmed by the stamp/ signature of the dealer on the warranty card.
- In cases of warranty or repair claims, please post the entire unit with the completed warranty card to the customer service point responsible for you.

- We reserve the right to make technical and optical changes as well as changes to the equipment!

For repair or maintenance, send the **Promed T-6** via your dealer to:

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The following information is vital for problem-free processing:

1. Original sales slip/receipt or dealer's stamp with the date of purchase
2. Defect detected
3. Unit name / type

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**Symbol for WARNING! Indicates sections of the user manual that are safety-relevant!**

The symbol corresponds to symbol 14, table D1, annex D of the EN 60601-1:1990 Standard. Instructions for a unit, system or accessories contain important safety-relevant information for the consumer.

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**Symbol for Type BF (body floating) units**

The symbol corresponds to symbol 2, table D2, annex D of the EN 60601-1:1990 Standard. Appliances that comply with the requirements of this norm offer a higher level of protection against electric shocks, especially when taking into account the permissible patient leakage currents and safety of the user.

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**Pay attention to the symbol for the user manual!**

The symbol corresponds to the D-M018 symbol in the DIN 4844-2 Standard. It is imperative that the manual be read before using the unit.

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**Symbol for Manufacturer and manufacturing Factory**

The symbol corresponds to symbol 5.12 of the EN 980:2008 Standard. It shows the name and the address of the manufacturer and the manufacturing factory.

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**Symbol for date of manufacture**

The symbol corresponds to symbol 5.6 of the EN 980:2008 Standard. It shows the date of manufacture.

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**Symbol for Notified Body**

This Symbol Complies with the European Medical Device Directive 93/42/EEC and amended by directive 2007/47/EEC requirements. Notified body TÜV Rhineland (0197).

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**Symbol for Authorized representative in the European Community**

The symbol corresponds to symbol 5.13 of the EN 980:2008 Standard. It shows the Name of Authorized representative in the European Community.

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**EN** | The requirements of the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EEC for class IIa medical products has been fulfilled according to the Standard EN 60601-1-2:2007.

CE according to the EC Directive 2004/108/EEC and EC Directives 2002/95/EEC, 1907/2006/EEC, 2005/69/EEC and 93/68/EEC.  
The CE marking for the appliance refers to the EC Directive 93/42/EEC in conjunction with EC Directive 2007/47/EEC.  
Appliance classification: Class IIa.

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Electromagnetic disturbances are neither emitted from the appliance nor affect the appliance.

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