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sanowell ComfortTherm Heizkissen

DE

EN

FR

● sanowell

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Thank you and congratulations.

We are happy that you have chosen a product by hofmann gmbh. Our name stands for innovative and high-quality products.

We want you to enjoy our products and to achieve the desired results. Please read the following instructions for use and care before using your electric heating pad.

The products of hofmann gmbh are subject to constant quality control. If you do find defects, or want to share suggestions or comments on this product, please do not hesitate to contact us.

You can reach us any time.

Your hofmann team

Delivered items:

The product package contains the following items:

- 1 sanowell ComfortTherm electric heating pad
- 1 cover for the sanowell ComfortTherm pad
- 1 Velcro fastening strip
- 1 manual in German, English and French

Safety instructions:

- Please read the manual carefully before using your heating pad for the first time.
- Do not use the heating pad while you sleep or if you are about to fall asleep, as prolonged use, or use at too high a setting, can lead to burns. Do not let children use the heating pad without supervision, keep the pad away from animals and never leave it unattended.
- Do not use your sanowell heating pad if you are sensitive to heat or if you are in a wet environment (e.g. the bathroom).
- Operate the heating pad only at 220~240V.
- Since the electromagnetic fields created by the product may lead to undesired interaction with pacemakers, consult your physician before using the pad. You should also ask a physician before using the heating pad on sensitive, irritated or inflamed skin.
- Check the heating pad for damage before every use. Such damage may be caused by careless handling, i.e. using the pad while it is bent or folded, damaging the pad with needles or other pointed objects, using the pad without the original cover, applying pressure to the pad (e.g. by sitting on it) and pulling or twisting the power cord.
- The switch may grow hot during operation and must therefore not be covered or placed on the heating pad.
- If you detect a defect in your sanowell heating pad, or if you have dropped the pad in water, pull the power plug out of the power outlet immediately. Do not touch the heating pad until the power circuit has been interrupted!
- Your sanowell heating pad may only be repaired and opened by an expert authorized by hofmann gmbh, otherwise the guarantee will expire.

Operation

With your sanowell electric heating pad, you can enjoy pleasant, relaxing warmth at several temperature settings. The fastening strip allows you to position the pad on the desired body part, applying heat and increasing the blood circulation in that area to soothe tense muscles.

Using your sanowell ComfortTherm S heating pad

- Plug the power plug into a power outlet.
- Switch on the heating pad by moving the temperature controller to the desired heat setting (1-4). The red lamp shows that the pad is heating up. You can regulate the temperature by moving the controller. 1 is the lowest, 4 the highest setting.
- After 15 minutes at setting 1 / 90 minutes at settings 2-4, the heating pad is shut off automatically. If you want to shut off the pad sooner, or if you want to start the pad again after the 15 / 90 minutes have elapsed, move the temperature controller to 0 and then after approx. 1 minute back to the desired heat setting.
- Switch off the pad by moving the controller back to "0". The red lamp goes off.
- Please remember to always remove the power plug from the power outlet after use.

Using your sanowell ComfortTherm L, XXL and Soft heating pad

- Plug the power plug into a power outlet. The orange lamp shows that the heating pad is ready for use.
- Use the On/Off button to switch on your sanowell heating pad. The "L" in the display shows that the pad is set to the lowest heat setting. To increase the heat, press the temperature button until the desired heat setting is shown in the display.
- After 180 minutes, the heating pad is shut off automatically. If you want to shut off the pad sooner, or if you want to start the pad again after the 180 minutes have elapsed, press the On/Off button.
- Please remember to always remove the power plug from the power outlet after use.

Cleaning and maintenance

- Pull the power plug and allow the pad to cool off, then clean it with a soft, slightly moist cloth or a soft brush. Make sure that no water penetrates the pad and do not use the pad until it is completely dry.
- Do not use aggressive cleaners. Store the pad in a cool, dry place; do not fold the pad.
- The cover, not the pad itself, can be washed as indicated by the symbols on the label.

Recycling

The Packaging should be disposed of in the garbage recycling, the metal parts of the device in the used metal collection. Plastic parts, electronic pieces and conductors need to be disposed of as electronic garbage. For further information concerning disposal of problem waste contact your nearest recycling center.

Specifications:

Model	ComfortTherm S	ComfortTherm L	ComfortTherm XXL	ComfortTherm Soft
Power supply:	220~240V, 50 Hz	220~240V, 50 Hz	220~240V, 50 Hz	220~240V, 50 Hz
Output:	60 watts	60 watts	60 watts	60 watts
Switch:	mechanical, 4 settings	digital, 11 settings	digital, 11 settings	digital, 11 settings
Automatic Shutoff:	Setting 1 after 15 minutes Setting 2- 4 after 90 minutes	after 180 minutes	after 180 minutes	after 180 minutes
Cord length:	approx. 200 cm	approx. 200 cm	approx. 200 cm	approx. 200 cm
Dimensions:	30 x 40 cm	35 x 50 cm	40 x 60 cm	35 x 45 cm
Weight:	0.45 kg	0.56 kg	0.63 kg	0.40 kg
Article-No:	89160	89161	89162	89163

Subject to technical and design changes in the course of product development.

Warranty

There is a 2 year guarantee for our products, applicable from date of purchase.

Within the prescribed warranty time we will remove all defects of the device, that are based on faulty material or manufacture. Removal of defect occurs at our choice by either repair or replacement of parts or the entire device.

Exempt from warranty are: Defects caused by inappropriate use (f. ex. connection to inappropriate electrical current or voltage, connection to inappropriate sources of electricity, breakage by fall and such). Normal wear and tear that only minimally affects the operational utility of the unit is not subject to the warranty.

Warranty subsides if unauthorized intervention or repair is undertaken, or if non-original parts are used.

In the case of a warranty claim please forward the product with the original receipt to one of our service partners or directly to the hofmann gmbh.

No other claims are valid under the provisions of this warranty, unless our liability is based on legal statute. Claims against the dealer based on the provision of the purchase agreement do not affect this warranty.

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