Gebrauchsanleitung | Instruction leaflet | Instruction d'utilisation | Instrucciones de empleo | Istruzioni per l'uso
Gebruiksaanwijzing | Руководство по эксплуатации | Instrukcja użytkowania | Bruksanvisning | Käyttöohje

promed 1020

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SE

FI

- A.** Socket for the power cord
- B.** Socket for the handpiece cable
- C.** Grinding head holder
- D.** ON/OFF switch
- E.** Speed controller
- F.** Rotation direction switch + standby
- G.** Holding device for the handpiece

Please position the station onto 0 or switch it off before changing the grinding head.

- 1.** Turn the knurl in an anticlockwise direction. The chuck opens.
- 2.** The grinding head can be changed.
- 3.** Following this, turn the knurl in a clockwise direction. The chuck closes and fixates the grinding head. Switch the unit on again for operation.

DESCRIPTION OF THE UNIT

Please keep this user manual in a safe place!

Please read the user manual thoroughly before using the appliance for the first time.

Thank you for purchasing the **Promed 1020** manicure and pedicure unit. This unit has been designed and manufactured in accordance with the Medical Device Directive 93/42/EEC and its amendments to guarantee quality of treatment. Both qualified specialists and private persons may use the unit after they have read the user manual. This unit is only suitable for use with manicures and gentle pedicures (e.g., removal of thickened skin and calluses, cuticle care, nail care, treatment of thickened nails, etc.). We, as the manufacturer, cannot be made liable in any way for injury or damages to people or objects that arise from failure to comply with this warning. The unit is also suitable for use in the area of nail techniques and nail design as well as the handling of all materials used in this area such as gel and acrylic.

We wish you a lot of enjoyment with your new **Promed 1020** unit. We would like to familiarize you with the unit in the following sections. Please read the user manual thoroughly before using the unit for the first time.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

INFORMATION ON THE UNIT / ACCESSORIES

Before operating the unit for the first time, please check that all contents have been delivered and that the line voltage is correct. This is done by checking whether the line voltage required for the unit (you can find this information on the specification plate) corresponds to the local power supply (operating voltage of 230 V ~ connected to a correctly installed wall socket with a supply voltage of 230 V ~; operating voltage of 115 V ~ connected to a correctly installed wall socket with a supply voltage of 115 V ~).

- The **Promed 1020** unit, in combination with the grinding heads designed for this unit, is ideal for manicures and gentle pedicures, and is also particularly suitable for use with nail techniques (nail design).
- The **Promed 1020** unit uses a high performance motor. The electronic speed control is used to smoothly set the speed of the grinding head. The unit can rotate clockwise or anticlockwise (reversible).
- The grinding head holder is able to fixate all commercial

grinding heads with standardized dental size shafts of 2.32 to 2.35 mm (0.091 - 0.093 inches). Please use only high-quality grinding heads with a rounded shaft end and a maximal shaft length of 32 mm (1.26 inches) to guarantee the longest life possible for the grinding head holder! Always make sure that the grinding head does not come in contact with other objects when you lay the handpiece down. There is a risk of the shaft becoming compressed and asymmetrical. Asymmetry in the grinding head causes damage to the handpiece!

- Grinding head bits are available at your local specialist retailer or **Promed**. (www.promed.de).
- In order to guarantee hygienic treatment, it is recommended that grinding heads be disinfected after use.
- We recommend that only wipe disinfectant be used for handpiece disinfection.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

OPERATING THE UNIT

Before you switch the unit on, you must set the rotating speed (rotary knob **(E)**) to zero.

The power plug is intended as a separator. Make sure that the power plug is plugged into an easily reachable wall socket.

- A.** Before operating this unit, please ensure that the voltage specified on the unit corresponds to the local line voltage. Plug the power plug into the wall socket and the other end of the power cord, found on the left-hand side of the unit, into its designated socket **(A)**.
- B.** Plug the handpiece's plug into the socket **(B)** on the right-hand side of the unit.
- C.** Position the rotation direction switch on "Standby" before using the grinding head. Select a suitable grinding head. Open the chuck **(C)** by rotating it counterclockwise. Insert the grinding head into its holder at the front of the handpiece. Lock and secure the grinding head by rotating the chuck

clockwise. Never use force to remove or insert a grinding head! Grinding heads with a shaft diameter of 2.32 mm to 2.35 mm (0.091 to 0.093 inches) fit exactly into the grinding head holder without you having to use excessive force. To remove the grinding head, please switch the unit off. Rotate the chuck in an anti-clockwise direction and remove the grinding head.

- D.** Push the toggle switch labeled "Power" **(D)** onto the "ON" position to switch the unit on. To switch the unit off, push the same switch onto the "OFF" position.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

OPERATING THE UNIT

E. Adjusting the speed:

The desired speed is set using the large rotary knob **(E)**.

F. Adjusting the direction of rotation:

To protect the motor and ball bearings, you should set the rotating speed to "0" or put the unit in "standby" for a short moment before you change the direction of rotation.

Clockwise:

To switch the unit into clockwise rotation, push the toggle switch labeled "Direction" **(F)** onto the "right" position.

Counterclockwise:

To switch the unit into counterclockwise rotation, push the toggle switch labeled "Direction" **(F)** onto the "left" position.

G. Standby mode: When the toggle switch **(F)** is positioned in the middle, the unit is in standby and the motor stops, e.g., to stop the grinding head for a short moment without altering the selected speed. However, to completely switch the unit off, please place the ON/OFF switch **(D)** onto the "OFF" position.

In cases of malfunctioning:

- Check that the power supply is correct.
- Check whether the handpiece is correctly connected to the unit.
- Check whether the unit is switched "ON".
- If the problems have still not been solved after following these steps, please consult your service center immediately.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

SPECIAL QUALITY FEATURES

ASC: Automatic speed control (digital load equalizing control):

The electronics register the pressure on the grinding head during grinding and balances out the loss of speed. You always have access to the maximal force.

OPC: Overload protection control:

If the grinding head is blocked or stopped for too long, the electronics register this and switch the unit onto a lower, safer performance level. However, you should not keep the unit in this state for a long period of time as this may damage the electronics.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

MAINTENANCE TIPS

- A.** Keep the inner parts of the handpiece as free as possible from dust and grinding residue.
- B.** Remove dirt from the outer surface of the handpiece with a cloth moistened with alcohol. Prevent moisture from entering the unit.
- C.** Only use grinding heads with a dental standard size of 2.32 - 2.35 mm (0.091 - 0.093 inches) and avoid using rusted or damaged grinding heads. Their false usage may lead to excessive wear and tear in the ball bearings. Adhere to the safe speed and the safety measures stipulated by the manufacturer.
- D.** Do not expose the handpiece to strong knocks as this may damage the ball bearings. If the handpiece should fall down, check whether the grinding head is bent. Pay closer attention to noise and/or heat generation. Send the handpiece immediately in for repair should it be damaged.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

Please pay attention to the following before first time operation:

Never place the unit in direct contact with fire, gas, oxygen or hot objects such as hot plates.

Please observe the following during operation:

- A.** The handpiece may only be connected to the unit included in the delivery contents.
- B.** Grinding heads to be used:
 - Standardized dental size: 2.32 - 2.35 mm / 0.091 - 0.093 inches.
 - Only use centered grinding heads with a rounded shaft end.
 - The use of rusted or bent grinding heads is not permitted.
- C.** Use the speed controller carefully.
- D.** Switch the unit off before you change the direction of rotation.
- E.** Before you switch the unit on, you must set the rotating speed (rotary knob **(E)**) to zero.
- F.** Appliance power supply: 115 V~ or 230 V~.
- G.** Switch the unit off before you change the grinding head.

- H.** Take every precaution possible to ensure that the unit does not fall down or get damaged in any other way.
- I.** If problems do arise with the handpiece, take it immediately to get it repaired.
- J.** Do not lubricate or wash the handpiece.
- K.** Do not switch the unit on while the chuck is opened.
- L.** Place the test grinding head in the chuck when storing the handpiece.
- M.** Never remove the dust covers from the switches as we can otherwise no longer guarantee that your unit will function properly.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

DANGER!

- Do not place the unit in direct contact with water or other fluids.
- Do not touch the unit with wet hands.
- Only operate the unit with the specified line voltage.

WARNING!

- Never leave the unit unattended when children or inexperienced persons are in its vicinity.
- Only use the unit for the applications described in this manual.
- Only use attachments and accessories from the manufacturer.
- Never use the unit after the electrical lead has been damaged, when the unit has fallen down or when it has come in contact with water.
- This unit is not designed to be used by people (including children) who have limited physical, sensory, or intellectual

abilities or who do not have sufficient experience and/or sufficient knowledge. This does not apply if these people are supervised by a person responsible for their safety or receive instructions from them on how to use this unit.

- Never let children play with the packaging material; there is a risk of suffocation.
- An electrical appliance is not a toy. Children do not recognize the risks that could develop when using electrical appliances. Therefore, operate and store the unit out of the reach of children.
- Do not place or use the unit on hot surfaces or near naked gas flames to ensure that the casing does not melt.
- Never let the unit come into contact with water. Use a moist cloth to clean the unit. Do not use solvents.
- Never carry or pull the unit by the power cord. Use the plug to pull the power cord out of the wall socket. Do not pull on the power cord itself.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

SAFETY INSTRUCTIONS

- Do not pull the power cord over sharp edges, wedge it in between objects, hang it down or position or hang it over hot objects (e.g., hot plates) or naked flames. Protect the power cord from heat or oil.
- Only use extension cords when you have satisfied yourself beforehand that they are not defect or damaged.
- If the unit's power cord should become damaged, always replace it with a new one. (Available in electrical specialist stores).

TRANSPORT/STORAGE/OPERATING CONDITIONS

- Only transport in the original packaging of the manufacturer.
- Transport and storage temperature:
+5° C – +50° C (+41° F – +122° F)
- Operating temperature:
+10° C – +40° C (+50° F – +104° F)
- Relative humidity: 30% to 75% non-condensing
- Atmospheric pressure: 700 to 1060 hPa
- An abrupt change in temperature is to be avoided as condensation may arise in the unit. Only operate the unit after the temperature has been equalized.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

SERVICE

Your customers will definitely be expecting you to deliver a professional result. This is only possible with an appliance that is neither defect nor damaged. Therefore, take advantage of our service and have your appliance regularly examined by our specialists. Send in your appliance to Customer Service every 12 months at the latest. We would be happy to provide you with a quote for the costs charged for appliance servicing.

DISPOSAL

Power tools, accessories, and packaging should be recycled in an environmentally suitable manner.

Only for EU countries:

Do not dispose of power tools with the domestic waste! In accordance with the European Directive 2002/96/EEC on waste electrical and electronic equipment and its implementation into national law, nonfunctioning power tools must be collected separately and recycled in an environmentally suitable manner.

DE

US

FR

IT

ES

NL

RU

PL

SE

FI

WARRANTY

This unit has been manufactured with the utmost of care and has been examined in detail before leaving the factory. We therefore provide a warranty of 24 months after the date of purchase when presented with the warranty card issued with this appliance, in accordance with the following conditions.

- Within the warranty period, we will replace all defect parts of the appliance free of charge, including labor costs for repairs under warranty, in cases of verifiable faults in material or manufacture that have occurred with correct usage and have been identified during the warranty period.
- This warranty is void when the unit has been improperly used (e.g., connection to unsuitable sources of electricity, breakage), the unit has been tampered with (e.g., opening the unit's casing), and when spare parts have been used that have not been approved of by **Promed**. Parts subject to wear are exempt from the warranty. Appliance motors and mobile parts are not included in the warranty.
- The warranty period begins on the date of purchase. Warranty claims do not affect the warranty duration. Warranty claims must be lodged within the warranty period. Claims arising after cessation of the warranty period cannot be considered.
- The warranty only comes into effect within the scope of these warranty conditions when the date of purchase has been confirmed by the stamp/ signature of the dealer on the warranty card.
- In cases of warranty or repair claims, please post the complete appliance with the completed warranty card to the customer service point responsible for you. The warranty card can be found at the end of this user manual.
- We reserve the right to make technical and optical changes as well as changes to the equipment!

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Symbol for protection class II appliances

The symbol corresponds to symbol 10, table D1, annex D of the EN 60601-1:1990 Standard. The appliance is double insulated as a protective measure. This prevents conductive parts of the appliance from becoming impaired, something that could lead to malfunctioning of the basic insulation voltage.

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Symbol for type BF (body floating) appliances

The symbol corresponds to symbol 2, table D2, annex D of the EN 60601-1:1990 Standard. Appliances that correspond to the requirements of this norm offer a higher level of protection against electric shocks, especially when taking into account the permissible patient leakage currents and the safety of the user.

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Symbol for **WARNING! Indicates sections of the user manual that are relevant for safety!**

The symbol corresponds to symbol 14, table D1, annex D of the EN 60601-1:1990 Standard. Instructions for an appliance, system or accessory that contain important information for the user that are relevant for safety.

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US | CE-Directives

The appliance is equipped with radio interference suppression in accordance with the EC Directive 87/308/EEC. It is designed for operation with a line voltage of 230/115 V~, 50/60 Hz. CE according to the EC Directive 2006/95/EEC and in accordance with the EC Directive 2004/108/EEC.

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The requirements of the EC Directive 93/42/EEC for class IIa medical products has been fulfilled according to the Standard EN 60601-1:1988 +A1:1991 +A2:1995; EN 60601-1:2003. The CE marking for the appliance refers to the EC Directive 93/42/EEC. Appliance classification: Class II. In accordance with EMC testing in EN 60601-1-2:1994; EN 55011 KLB:1998; EN 61000-4-2:1996; EN 61000-4-3:1997; electromagnetic disturbances are neither emitted from the appliance nor affect the appliance.

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Promed GmbH